

TA[®]iQ Rear Wheel Drive

Enable M.O.H
Panel Approved

When you require optimum traction



Outdoor
Performance



Sit under any Table



Reach Higher

Navigates easily around the home and through shops. Sits lower under standard dining tables.

Rides with ease through parks, streets and over curbs.

Elevates to give increased reach to benches, appliances, storage or to socialise.

All TA chairs come standard with Tilt and a HiLo function



Morton & Perry
Healthcare Equipment Solutions



RWD



HiLo function comes as standard with the TA RWD 30cm/12" elevation



Extra Low only 38cm/15" from floor to seat base

Benefits versus other powerchairs

Lower seat height	Sit under any standard table or desk
HiLo function	Reach benches and cupboards
Optimum Range	40km on full charge
Safety in transport	Crash test approved safe for user to remain in chair

RWD Outdoor Performance

The majority of the weight in the RWD is centered around the drive wheels at the rear. With only 4 wheels this provides improved power and torque ideal in more difficult terrain.

The lowest standard floor to seat height on the market

A very low floor-to-seat height (excluding cushion) of 38cm/15", achieved with HiLo function, allows the user to sit at tables and desks more easily. The low seat height also increases visibility and head clearance when travelling in a vehicle.

Higher electric elevation

Higher elevation of 30cm/12" (38-70cm/15-27½") ensures basins and higher cupboards are within easier reach.

High Performance

With a maximum programmable speed of 12.5km/h, the user can travel quicker when they need to, over roads or to avoid the rain. The TA iQ has an optimum range of 40km.

Can be configured easily

For all user heights from 3'6"-6'6"/107-198cm and related leg length without compromising ground clearance and speed.

Model	TAiQ RWD
Drive configuration	RWD
Product code	XTAIQRWD02
SAFETY Compliant - to ISO7176-19, NZ LVVTA 45-60 2.6 Wheelchair design	
Occupant crash tested	Yes
FUNCTION *Standard specification, no up-charge	
Seat tilt, electric	45° *
Seat lift, electric	30cm/12" *
Floor to seatpan minimum	38-43cm/15-17"
Floor to seatpan maximum	68-73cm/27-29"
Back recline, manual or	50° *
Back recline, electric	50° *
Legrest centre, manual	Yes *
Legrest swing-away, manual	Optional *
Legrest elevating, manual	Optional
Legrest centre, electric	Optional
Legrest swing-away, electric	Optional
Seat glide - 10cm forward	Optional
Specialty control options	Optional
PERFORMANCE	
Turning radius	65cm/25.5"
Chair width	60-63cm/23.5-25"
Chair length without footrests	91.5cm/36"
Range max.	40km
Speed max.	12.5km/h
Client weight max.	140Kg
SPECIFICATIONS	
Motor type	4 Pole
Battery size	80Ah x 2
Battery charger Amp	8Amp
Electronics, Penny & Giles	RNET expandable
Drive wheel size	35cm/14"
Castor front/rear size	20cm/8"
Powerchair weight	140Kg
ADJUSTABILITY	
Seat width	37-55cm/14.5-21.5"
Backrest width	37-55cm/14.5-21.5"
Effective seat depth	35-55cm/14-22"



Morton & Perry
Healthcare Equipment Solutions

Contact our team for a product demonstration
0800 238 523 | mortonperry.co.nz