




THE WORLD'S FIRST

Total shear measurement device

ishear.com

THE WORLD'S FIRST

The first clinical tool to measure
total shear force in the seat plane



ORDER
NOW
From
Morton & Perry

Setup is easy

All you need is iShear
and a mobile device.
Simply download the
free iShear app and
off you go.
It's shear genius.



SHEAR GENIUS.

iShear helps you make seating outstanding

Sliding is a huge clinical problem but difficult to quantify. Until now! Introducing iShear: the smart device that measures Total Shear Force*, helping you evaluate the positioning of a person in a seating system.

The ultimate solution for set-up and long-term use

iShear helps you minimize sliding forces when setting up the wheelchair. iShear also allows you to educate wheelchair users and care givers about the ideal use of the wheelchair.



* **Total Shear Force (TSF):** The total shear force in the seat plane measured by the iShear is the anterior-posterior force parallel to the seat plane. It's measured between the cushion and the seat base interface.

Put your clients in a great position



iShear device

- For wheelchairs from 35 up to 50 cm (width and depth)
- Bluetooth connection between iShear and mobile device
- Up to 8h over-time logging
- iShear is placed underneath the cushion

iShear App

- Real-time measurements
- Save, export and share measurements
- Compare multiple measurements
- Add notes per measurement
- Free download App



Real-time view



Show files



Compare



Log session

DOWNLOAD
FREE APP



Morton & Perry
Healthcare Equipment Solutions