

# Alternating Mattress



BetterLiving alternating pressure surfaces are designed to offer a cost-effective solution for the prevention and management of pressure related injuries.

Harnessing the latest technical features for improved clinical care and ease of use, BetterLiving offers a mattress replacement ideally suited to home care or aged care environments, this products can be used for people at risk to high risk of pressure injury.

<b>ITEM</b>	Mattress Replacement
<b>CODE</b>	APMBL-R01
<b>FEATURES</b>	<ul style="list-style-type: none"><li>○ Essential features for pressure injury prevention and treatment</li><li>○ Automatic pressure setting</li><li>○ Replace existing mattress for superior therapy</li><li>○ Cell in cell for enhanced user safety and comfort</li><li>○ Easy to operate and maintain</li></ul>
<b>CAPACITY</b>	40 - 180 kg
<b>PRESSURE INJURY RISK</b>	High to very high risk
<b>COMPLIANCE</b>	IEC 60601-1, IEC60601-1-2, IEC60601-1-11
<b>ARTG</b>	289458
<b>WARRANTY</b>	1 year

## Why choose BetterLiving?



### Safety

No more concerns over correct pressure level settings with our AutoCair feature – automatic weight sensing technology that delivers fail safe patient comfort levels irrespective of their position.

With the AutoCair feature active, the system adjusts to the appropriate pressure setting based on the patient's weight, reducing the need for carer intervention. Manual adjustments can still be made for individualised therapy or comfort.

### Simplicity

Offering all the features you need, the BetterLiving pressure surfaces are designed for basic ease of use and minimal maintenance. Easy to set up, simple to operate and straightforward to clean and maintain, these products require minimal technical training or support.

### Quality

Our BetterLiving pressure surfaces conform to all relevant Australian standards. Made using the best quality materials, every product must pass our stringent Quality Control inspection and regular testing protocols.

