



Workplace Chairs



4 Reasons to consider

1. Pain and posture management
2. Improved sitting endurance
3. Mobility and stability in the seated position
4. Adjustable for different working heights



Morton & Perry
Healthcare Equipment Solutions

0800 238 523 | mortonperry.co.nz



Benefits over a regular office chair

Greater Stability & Mobility

- > The sturdy wheel brake keeps the chair anchored so it will not suddenly move.
- > A stable chair can be relied upon to assist when lowering to sit or pushing up into standing.
- > The four wide & soft roll wheels allow for straight line foot propulsion or easy glide between desk and surrounding surfaces using your arms. Reducing the need to stand.
- > The stable base provides safe sideways active reach within the working space without the danger of tipping.
- > Move the chair with ease into or away from the workstation, swivel the base using your feet or release the seat lever to rotate it 90 degrees from the work station.

Adjustable height

- > Transfer in or out of the chair with more confidence & safety.
- > The electric lift option can save energy with standing from sitting.
- > Built from high quality materials the chair has a power or manual adjustable height function.
- > Combine the change in height with altering the tilt on the seat-pan to keep your feet in best contact with the floor even if you use different work stations, sitting or perching.
- > Use the brake, extensive range of height adjustment and seat tilt to reduce the work of standing. By starting the movement from lower height with your feet on the floor and body in a ready position.
- > Reach higher from sitting if needed.

Improve Posture - Reduce Pain

- > Vela chairs have a range of seat angle and back angle adjustments.
- > Having postural variety can assist pain and posture by allowing many varied sitting positions.
- > Powered elevate can increase weight bearing through the feet without having to actually stand.
- > Powered elevate can reduce the stress and strain of injured tissues by taking the work out of sit to stand .
- > Being able to easily move the chair within the work station by feet or arm propulsion can reduce the number of sit to stand actions needed at work.

Increased Endurance

- > Improved endurance for sitting and working can be achieved.
- > The Vela chair is safe and easy to sit into and get out of.
- > Postural variety via seat functions allows you to keep sitting but change the position of your body to ease overuse
- > Safe reaching from a stable chair can reduce twisting
- > Easy mobility allows more efficient work.

Electric HiLo Function



Sit to stand with electric HiLo adjustment with a simple push of the button

Walkable



Around the office.

Ideal for the following user conditions:

Health

- > Back pain

Neurological

- > Balance problems
- > Parkinsonism

Musculoskeletal

- > Muscular dystrophy
- > MND
- > Multiple sclerosis

Trauma & Medical

- > Hemiplegia
- > Spinal disorder
- > Cerebral palsy
- > Amputation
- > Restricted growth

- > Mobility impairment

- > Walking impairment

- > Arthritis

- > Stroke

- > COPD (Lung Disease)

Tango 700E

Close-up



→ **Backrest**

- > 'Active' or 'Contour'
- > Tilt: 10° Forward, 15° Backward

→ **Armrest**

- > Width & height adjustable
- > Foldable armrest optional

Brake Lever ←

Right hand or Left hand

→ **Seat**

- > 'Active' or 'Contour'
- > Tilt: 10° Forward, 15° Backward
- > Rotatable 90° to the left or right

Height Adjustment

Electric 20cm

Castors

Four big smooth rolling 10cm high x 4cm wide castors

Seating Customisation options

We will build each chair to the user's unique fit and function requirements.



Active

Offers great comfort and light support. It comes in various sizes to fit the user's body and support needs.

Seat available in: Small, Medium, Large, XLarge & Deep

Backrest available in: Small, Medium, Large & XLarge



Contour

Offers excellent comfort and a high degree of support. The seat is shaped with side support that helps stabilise your sitting position.

Seat available in: Small, Medium, Large & XLarge

Backrest available in: Medium, Large & XLarge



Foldable armrest



Brake extender with T-Lever



Anti tilt bracket



Footrest fold-up with legs



Footrest Swing-away set

Electrical options (INDENT): Sit-to-Stand Assist and Electric Brake.
If you require other customisations please **Talk to us for more options.**

Work safer with improved energy levels

The brake keeps you stable while sitting. A regular office chair can't guarantee this, and it can be exhausting to waste your strength and energy worrying whether the chair is going to roll away.

You can 'walk' or push a VELA mobility chair around while sitting down, and that means you don't have to waste unnecessary energy on getting up and sitting down – you just 'walk'.

If you need to grab something in the top cupboard at your workplace, the chair's electric lift can bring you up to the right height. The brake ensures that the chair doesn't roll out from under you while raising or lowering the chair's height – it's safe and secure.

Furthermore, the chair's comfortable backrest and seat, as well as the chair's ability to be adjusted to fit you ergonomically, can help reduce your pain.



If you sit at a desk every day, a work chair from VELA can provide excellent conditions for changing your working posture during the day.

The seat can be angled forward so that your feet have a good support surface, and you get an open hip angle that makes it easier to sit nicely at a desk. In addition, the proper sitting height and sitting posture mean that your shoulders automatically become more relaxed, providing a good working posture.

For raised desks lock the chair's wheels and tilt the seat forward. It supports a partially standing working position, enabling you to have varied working postures during a long day. The chair still supports you, and it's not too stressful for your legs and back.



Leanne from Southland

Medical Condition: Muscular Dystrophy

Primary Use: Office

The Vela Chair Allowed Leanne to Continue to Work

"I have been using this chair in an office setting since early 2020. The difference this chair has made to me is between continuing to work or having to leave my job. That to me is a big deal, as it is hard enough to find employment when you have a disability. Prior to having this chair my legs and feet were extremely tired and sore by the end of the day and then very restless at night making it hard to sleep. The ability for the chair to lower right down to nearly the floor, has been very helpful to get to the bottom draws of the filing cabinets etc. I no longer have to ask someone for help all the time. Being as independent as possible is a big plus."

Exclusively in New Zealand available from:



Contact us for a product demonstration:

0800 238 523 | mortonperry.co.nz