


# ProCair Prime

The ground-breaking, unique alternating pressure mattress with an integrated pump.




Base Cover with dual zip for easy removal




Integrated pump

Integrated pump eliminates costly pump repairs from pumps being damaged when bumped onto the floor



Can be adapted between dynamic mode or hybrid mode based on different clinical needs



Automatic weight sensing technology delivers optimal pressure care

Easy to setup, the ProCair Prime automatically adjusts the pressure to suit the weight of the patient and does not need a separate pump unit to hang from the bed - it is all built in!



The ProCair Prime is adaptable for different clinical needs. It can be used as a fully alternating pressure mattress (dynamic mode) or easily converted to hybrid mode by positioning the foam sheet above the air cell layer.

## SPECIFICATIONS

<b>CODE</b>	XNAPMPC-R03
<b>CAPACITY (SAFE WORKING LOAD)</b>	30-200 kg
<b>SYSTEM ALTERNATION</b>	1 in 2
<b>WARRANTY</b>	3 year
<b>PRESSURE INJURY RISK</b>	Very high risk
<b>COMPLIANCE</b>	IEC 60601-1, IEC60601-1-2, IEC60601-1-11
<b>ARTG</b>	289458
<b>MATTRESS DIMENSIONS (L x W x H)</b>	2000 x 900 x 180mm
<b>MATTRESS WEIGHT</b>	12 kg



Refer to the Technical Data Sheet for detailed specifications.

PCPMB/0519 MAY 2019